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Service commercial :

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CONTENTS	2
I - INTRODUCTION	3
II - COMPONENTS	3
III - SETTLEMENT OF THE SENSOR AND OF THE BOX	4
IV - MUSCULAR SWITCH USE	4
V - PROBLEMS	5

CAUTION

THE 7N93 SWITCH IS NOT AN ALARM SYSTEM.

ITS USE PRINCIPALLY CONCERNS THE ELECTRONICAL TECHNICAL AIDS
CONTROL
(COMMUNICATION AIDS, ENVIRONMENTAL CONTROL, COMPUTER ACCESS,...)

DO NOT CONNECT ANY EQUIPMENT WHICH IS MORE THAN 24 V AND 0.25 AMPS

CONSIDERING THE SUPPLIED ELASTIC WEBBING CONTAINING RUBBER,
PLEASE VERIFY THAT IT DOES NOT PROVIDE YOU ANY ALLERGY ON YOUR SKIN.

THIS EQUIPMENT IS SUPPLIED WITH A BATTERY. PLEASE READ THE NOTE
CONCERNING THE BATTERIES.

IT IS ADVISED TO REMOVE THE BATTERY BEFORE TRANSPORTATION
OR IF YOU DO NOT USE IT DURING A LONG TIME.

PLEASE DO NOT WALK ON THE SENSOR.
DO NOT BEND THE SENSOR, WHICH WOULD HARM THE GOOD WORKING OF THE
SWITCH.

I - INTRODUCTION

The switch 7N93 is a polyvalent one, mounted with a sensor with little dimensions ; these little dimensions allow it to put it easily on the body. It is delivered with a Velcro headband which facilitates the installation on the forehead or on the arm for exemple.

There are two possible adjustments :

- sensitivity,
- the time between two actions.

The switch 7N93 detects the muscular contractions, but weak changes of temperature too.

II - COMPONENTS

- A control box with a label showing the serial number, and with different switches with two possible adjustments.
- A sensor which looks like a flat circle with a diameter of 1,5 cm for which the two faces are recover with Velcro. It is connected to a cable about 1 m long.
- A contactor cable which allows to connect the control box to an electrical equipment.
- A strap which can be used as a headband, bracelet, etc. to settle the sensor.

III - SETTLEMENT OF THE SENSOR AND OF THE BOX


Diagram 1 :
Upper face of the control box

Diagram 2 :
Front face of the control box

- Connect the sensor (the cable with its circle) on the sensor plug marked 
- Place the sensor on the skin (for exemple on the forehead with the headband supplied).
Caution ! The choice of the sensor place is very important : the movement has to be controlled and has to be regular, without any interference. The sensor can be fixed either with the supplied strap, or with the adhesive medical ribbon. Verify that the sensor is not too squeezed on the skin and that the sensor wire does not move, can not be easily drawn, which could open a contact.
- Place the interrupter on position 
- Adjust the sensitivity as well as the sensor react only to a voluntary movement (the switch makes a noise).
- Then, you can do a test in placing, for exemple, the sensor at the extremity of a finger with a medical adhesive rubbon, then in pressing on any surface with the finger. This manipulation converts the muscular switch in sensitive switch and gets it in permanent contact with the skin without any risk to loose this contact.
- As well as the test finished, it is possible to connect the muscular switch on an equipment to be controlled with the supplied cable.
- Connect this cable on the plug marked 

CAUTION : DO NOT EXCEED 24 VOLTS AND 0.25 AMPS
FOR THE CONNECTION OF AN EQUIPMENT TO THE MUSCULAR SWITCH

IV - MUSCULAR SWITCH USE

- Light diode of function :

The green light diode comes on during the functioning of the switch : this functioning is an impulse which is about 1 second. The muscular switch does not work continuously.

- Temporizing key marked

This key allows to choice if you wish a temporizing between each impulse or not. It can be useful in the case of parasite movements after the voluntary movement.

- Adjustment screw temporizing marked

This screw allows to personalize the temporizing from 0.5 sec to 2 sec. During all the temporizing, the red light diode remains switched on, indicating that an action is not possible.

- Interrupter Run/Stop

This interrupter has 3 positions :
a position stop at the middle
a position run with a bip
and a position run without bip.

V - PROBLEMS

The control box does not work, the green diode does not switch at the movement.

Cause	Solution
The battery is exhausted.....	Remote and change the battery.
The sensor is damaged.....	Verify the sensor and change it if necessary.
The cables are not well connected	Connect the cables correctly and verify the connectors.
The interrupter RUN/STOP has STOP position	Set RUN position.

The contactor works but it is not so sensitive

The battery is exhausted	Change the battery.
The adjustment of the sensitivity is not good ...	Modify the adjustment of the sensitivity.
The sensor is not well settled or is damaged ...	Settle it again or change it if necessary.

The switch works but it is too sensitive

The battery is exhausted	Change the battery.
The adjustment of the sensitivity is not good	Modify the adjustment of the sensitivity.

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