**ORTHÈSE MAG 2**

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**MAG 2 ORTHOSIS**

Reference 7A2099 - Page 16/29

**ORTESIS MAG 2**

Referencia 7A2099 - Page 30/43

**MAG 2 ORTHÈSE**

Anleitungsnummer 7A2099 - Page 44/57

Service commercial / Customers service / Servicio comercial /  
Kundenabteilung :

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## CONTENTS

---

<b>I - PRINCIPLES</b> .....	<b>17</b>
<b>II - MAG 2 SUMMARY CHART</b> .....	<b>18</b>
<b>III - COMPONENTS DESCRIPTION</b> .....	<b>19</b>
<b>IV - INSTALLATION - SETTINGS</b>	
A. Installation .....	<b>23</b>
B. Settings	
1. Pitch and roll .....	<b>23</b>
2. Range-of-motion limits .....	<b>24</b>
3. Yoke and adjustable balance point .....	<b>24</b>
<b>V - PREPARATORY WORK</b>	
A. Muscular, articular, sensory and proprioceptive examinations .....	<b>25</b>
B. Settings	
1. Pitch and roll angles .....	<b>25</b>
2. Rotation stop and slide stop locations .....	<b>26</b>
3. Adjustable yoke and balance point .....	<b>26</b>
C. Exercises and training .....	<b>26</b>
D. Specific case of utilization with electric wheelchair .....	<b>26</b>
<b>VI - APPLICATIONS</b> .....	<b>27</b>
<b>VII - APPENDIX</b>	
Wheelchair installation kit 7A09 for MAG 2 .....	<b>28</b>

*Thank you for buying the MAG 2 arm support orthosis.*

*MAG 2 is the product of research in close collaboration with the Daniel Medical-Teaching Centre of St-Hilaire du Touvet (France).*

*To get the most from your MAG 2 we recommend that you read this manual carefully.*

*You will see that MAG 2 can be adjusted in many ways and can be adapted to any situation.*

## **I - PRINCIPLES**

---

MAG 2 is an arm support orthosis.


It reduces “downpull” substantially and eliminates friction.

It enables you to use muscles that would otherwise be useless.

Shoulder movements are relayed to the hand by a system of levers and pivots.

It is entirely modular.

## MAG 2 SUMMARY CHART

*Nota : The MAG 2 orthosis - product number 7A20 - includes product numbers 7A21, 7A22, 7A23, 7A24 and 7A26. The other product numbers are for accessories.*


### III - COMPONENTS DESCRIPTION

---

#### Reference 7 A 21 - Tabletop attachment

This attachment is used to fit MAG 2 to a tabletop from 18 to 50 mm thick and to adjust its height.


*Figure 1*

#### Reference 7 A 22 - Tilt blocks


These two blocks rotate separately around a centre shaft. By rotating them, you can select the pitch and roll of your MAG 2.

*Figure 2*

### Reference 7 A 24 - Slide / Carriage Assembly


The original feature about your MAG 2 is the slide principle that allows both precise movement and a large range-of-motion. Being compact, it is unobtrusive in the rest position.

*Figure 3*

### Reference 7 A 25 - Short Connecting Rod (centre distance 70 mm)

The connecting rod fits between the tilt blocks and the slide-carriage assembly to provide a large range-of-motion for the arm.

*Figure 4*


## Reference 7 A 26 - Standard Cast / Adjustable Yoke Assembly

The standard small cast comes with velcro straps to maintain the forearm in position.

The cast is not designed for all uses. Optimal use of the MAG 2 requires the design of a customized cast by an ergotherapist or orthopedist.


The adjustable yoke restricts or completely prevents motion of the pivot. The pivot angle can be adjusted from 0 to  $\pm 45^\circ$ .

This yoke allows the point of equilibrium to be adjusted for movement with minimum effort.

*Figure 5*

**Reference 7 A 27 - Long Connecting Rod (centre distance 120 mm) : (optional)**

*Figure 6*

**Reference 7 A 09 - Wheelchair Installation Kit : (optional)**

This kit fits any diameter 15 - 25 mm round or square tube.

It can be adjusted in all three dimensions (length, width and height).

*See Appendix page 29 : Wheelchair Installation Kit.*

## IV - INSTALLATION / SETTINGS

---

### A. Installation

1. Your MAG 2 can be installed as a permanent unit on a tabletop or mounted on a wheelchair.

tabletop : use attachment 7 A 21, to which the tilt block 7 A 22 is fitted.

N.B. : To prevent premature wear, do not overtight the different adjustment screws.

#### Wheelchair

See Appendix page 29 : Wheelchair Installation Kit.

After fitting the installation kit onto the wheelchair, screw the two tilt blocks onto the square horizontal tube.

☞ It may prove important, after installing the MAG 2 on the wheelchair, to saw off the unused lengths of the two round and square tubes of the kit so they take up less space.

2. After conducting enough trials to determine the required range-of-motion for arm movements and part clearances, place the tilt blocks :
  - either directly on the slides 7 A 24,
  - or the short connecting rod 7 A 23,
  - or the long connecting rod 7 A 26,
  - or possibly the two connecting rods (short, then long or the opposite), followed by the slides.
3. Next fit the standard cast (7A 26) or the special cast on the carriage via the yoke.

### B. Settings

#### 1. Pitch and roll (Figure 2)

Your MAG 2 comes with a hex wrench. This wrench can be used to unscrew the tilt block retaining screws.

The blocks can then be rotated relative to one another to adjust the angle of pitch and rotated around their shaft to select the exact orientation required.

☞ This adjustment is particularly important for the future use of your MAG 2.  
*See page 26 : Medical Examinations.*

## 2. Range-of-motion limits (Figure 7)

The two connecting rods come with metal studs that can limit the rotation of the connecting rods and slide.

You will have noticed that the connecting rods and tilt blocks have a circle of 16 holes drilled in them. These holes accommodate the metal stops provided. They can limit rotation from 0° to 180° in 22.5° steps.

☞ The rotational amplitude setting prevents excess rotational motion from which it is impossible to return to the rest position.

See page 26 : B. Settings.


Figure 7

## 3. Yoke and adjustable balance point (Figure 5)

The Yoke is fitted with a pivot with two screws. The position of these screws adjusts the movement of the cast pivot very precisely.

☞ It is important to be able to adjust this movement, particularly during initial use. Users sometimes have difficulty in controlling the free pivot movement of the cast. Being able to adjust this motion very precisely makes it easier to adapt gradually.

The bar fixed on the underside of the cast can modify the balance point of the cast.

This accessory is very useful if your MAG 2 requires a customized cast. In this case, the balance point is very often modified. A well-balanced cast is easier to move.

## V - PREPARATORY WORK

---

### A. Muscular, articular, sensory and proprioceptive examinations

Before using your MAG 2 it is important to have :

- a muscular examination

This examination will indicate the results that can be expected from your MAG 2. Residual muscular activity in the shoulder (the deltoid in particular) is an important, but not necessary, starting point. The possibility of using the biceps, even if weak, is an asset.

- an articular examination

Flexible joints with a good range-of-motion, without locking and with no incorrect attitudes will make it easier to use the MAG 2.

- a sensory and in particular proprioceptive examination

For the best functional results.

A series of exercises in a rehabilitation centre, with an upper limb suspension appliance will make it easier to assess the potential for movement of the shoulder and elbow.

These exercises will also indicate the different compensation mechanisms to achieve optimum functional results.

### B. Settings

Once the MAG 2 system is installed, the next job is to define :

1 - the pitch angle and block orientation (*See page 24 : Installation-Settings*)

The choice of these two angles is very important. They must be chosen so that gravity takes the orthosis in the direction the user cannot go alone. As a rule, these angles are calculated so that the arm comes back easily to the neutral position close to the body.

Graduations on the tilt blocks 7 A 22 are used to locate the specific settings for each patient.

2 - the position of the rotation stops and slide stops

It is also important to restrict the range-of-motion. This limit must prevent the user from reaching a position from which recovery is impossible (e.g. excessive elbow extension or lateral rotation of the shoulder).

### 3 - the yoke and the balance point setting

So as to reduce as far as possible the effort required to keep the cast in position and to seek a balanced position.

Example :

Where a keyboard is used, the horizontal position of the cast must be correct. If muscle strength is inadequate, a low amplitude setting of the yoke pivot is required.

## **C. Exercises and training**

Two valuable exercises are advocated by the Daniel Douady Medical-Teaching Centre, the MAG 2 orthosis designers.

- Test on computer joystick with real tests on the electric wheelchair control.
- Drawing and graphic exercises on a worktop to determine the maximum functional amplitude. To carry out these exercises :
  - Set up the one-side or two-side MAG 2 orthosis.
  - Fasten a large sheet of paper to the tabletop.
  - Fix a pen or pencil, etc. to the cast.
  - Ask the user to draw the largest possible circles. These can be used to determine the zones reached by the hand or tool and to adapt them to the constraints of the activity.

## **D. Specific case of utilization with electric wheelchair**

*See Appendix page 29 : Wheelchair Installation Kit.*

- When adjusting the wheelchair installation kit, the position of the proximal shaft of the MAG 2 orthosis is important and special care must be taken with it. The shaft may be close to the hand or close to the elbow.
- For manipulating the wheelchair control joystick, if the standard cast is unsuitable, it is important that the customized cast be made so as to “ hold ” the wheelchair joystick (e.g. contracture of triceps in extension), but also to move away easily if needed.

One easy and common technique is to stick a non-slip film under the cast “palm”.

## VI - APPLICATIONS

---

There are many very positive applications in functional and psychological terms :

### **Functional benefits :**

- give amplitude to hand movements
- use the electric wheelchair standard joystick
- use a keyboard (computer, typewriter)
- turn the pages of a book
- feeding
- movement guide in the case of trembling or slight uncontrolled movements.

### **Rehabilitation :**

- strengthening of muscles that are not fully paralysed by multiple use every day, with sometimes one point gain in a matter of months
- auto-mobilization of the shoulder and scapular girdle with reduced contracture and possible disappearance of chronic scapular pain.

### **Psychological benefits :**

- Recovered mobility of the upper limbs changes the image the disabled person projects and restores expressive capability to the hands.

*Other accessories can be fitted. Please ask for further information.*

## RESTRICTIONS OF USE

---

The average weight of an arm rested on a support is about 3 kg, depending on the person.

The safety load of MAG 2 enabling a normal movement is 6 kg.

☞ MAG 2 was tested with a 4-kg weight placed on the standard cast, with all joints in motion. MAG 2 was tested in this configuration during 93 400 cycles.


☞ MAG 2 is mainly advisable for pathologies such as traumatic or neurologic tetraplegia, muscular dystrophy... It is not designed to guide movements in situations of ataxia.

## VII - APPENDIX : WHEELCHAIR INSTALLATION KIT 7A09 FOR MAG 2


This kit comprises :

- 2 shims with clamp screw,
- 1 obliquely drilled attachment,
- 1 round tube with T-piece to accommodate a square tube ; the tilt blocks fit onto the square tube.

1. Undo the 2 hex socket screws on each shim and fit them around the wheelchair tube. Retighten the two hex socket screws.


2. Insert the drilled adapter ① through each shim between the wheelchair tube and the clamp screw and then tighten it (see below).

3. Fit the round tube into the hole in adapter ① on the wheelchair.
4. Fit the square tube and fasten the two tilt blocks 7 A 03 of the MAG 2 onto the square tube.

Tilt blocks

The kit can be adjusted in all dimensions (length, width and height).

It is sometimes necessary to remove the wheelchair arm-rest or the arm-rest padding.

The 2 shims are used to clamp any round or square tube of 15-25 mm diameter or cross-section.

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